




Consider Them Your Personal Support Group



MST Capsule Retractors "Classic"



MST Capsule Retractors, Chang Modification

When we looked at and talked to people about what they liked and disliked about capsule retractors, we kept hearing the same thing, "I like the support they give, but the amount of focal pressure on the capsulorhexis margin is unsettling".

So... we set out to design a new kind of capsule retractor: one that puts less focal pressure on the capsulorhexis margin. We now provide you two options. 1. MST Capsule Retractors "classic" that distributes force in the equator and 2. A modified design by David Chang, MD that reduces the rare chance of a CTR being threaded in the open end of our classic capsule retractor.

MST Capsule Retractors. Only from MST.