

Catalog Number and Packaging

Provided sterile. 6 packages of 4 MST Capsule Retractors per box.

Description:
MST Capsule Retractors
MST Capsule Retractors,
Chang Modification


Catalog Number
MCR-0001
MCR-0002



There are 4 retractors in each package, six packages in a box.




The retractors are held in place by a silicone bar and are gripped by slots in the bar.




The Tyvek lid is pulled back and the tray turned upside down. The bar will drop out.



Like this...



The bar is picked up and squeezed as shown to release its grip on the retractors.



Remove the retractors.

