

REFRACTION UNIT TRU-800

SMART REFRACTION UNIT

DELIGHT IN SIGHT

Mini. Smart. For everybody.



- Mini + smart
- Easy to install
- Accessory shelf
- Projector column with overhead light
- Electronically controlled patient chair
- Adjustable phoropter arm

 **TOMEY**
TECHNOLOGY AND VISION

THE TOMEY REFRACTION UNIT TRU- 800



QUALITY IN DETAIL

Perfect if you need a tiny and clever solution for your refraction environment. Much smaller than standard refraction units does it even fit to very small rooms. For a proper refraction just attach all essential instruments like manual/automatic phoropter and chart projector to the TRU-800 . The control unit fits to the little shelf. If you want space for hand-held units or trial lenses - the optional trial lens cabinet TC-1000 meets with the design.



Accessory shelf with controller keyboard



TRU-800 foot switch (optional)



TRU-800 with standard chair SC-1000

To be combined with electronically reclinable chair or standard chair

SPECIFICATIONS

Shelf	Accessory shelf
Shelf height	75 cm
Seat height	44 cm (ER-1000 / SC-1000)
Seat elevation	44 – 59 cm (ER-1000 / SC-1000)
Seat reclination	30° (ER-1000)
Nominal lifting capacity	max 120 kg
Projector column height	170 cm
Overhead light	white LED 3.5 V
Power consumption	300 VA
Power supply	230 V
Weight	approx. 108 kg (ER-1000) approx. 85 kg (SC-1000)
Optionals: Phoropter arm, trial lens cabinet 61 x 61 x 79 cm (WxDxH), ophthalmoscope holder , foot switch	



Cabinet TC-1000