

# OASIS®



## *Corneal Comfort & Protection.*

- ▶ Ocular surface protection following surgery, injection, traumatic and non-traumatic corneal conditions
- ▶ Comfort of a thin, clear, and pliable collagen barrier
- ▶ Hydrates to 70% water content
- ▶ Protects the surface of the eye
- ▶ Lubricates the ocular surface with collagen while it gradually dissolves

Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times

Diameter: 14mm      Compound base curve: 9mm

## More than expected

OASIS® is dedicated to equipping eye care providers with effective solutions to ensure quality patient care, improve practice management, and support practice growth.

Find out more about how we can help your practice. Contact OASIS® for more information.

[www.oasismedical.com](http://www.oasismedical.com)

800-528-9786 (USA Toll Free)

909-305-5400

SOFT SHIELD and OASIS names and logos are registered trademarks of OASIS Medical, Inc.  
514 S. Vermont Ave, Glendora, CA 91741. LIT-SFSHC-S.Rev1-4.2014

