

surgiLine

The modern surgeon's stool.



OPHTHALMIC
SURGERY

Order No. 22000008

ADVANTAGES



Armrests

- Rugged, 360-degree movable, and height-adjustable
- Designed for personalized adjustment to the surgeon

Seat System


- Seat system evenly distributes your weight as you make the slightest adjustment in position.
- Anatomically formed, ergonomic backrest reduces pressure on inter vertebral discs and provides ideal support for the lumbar spine.
- Seat system upholstery in standard color anthracite F6461197
- Outstanding sitting comfort

Chassis

- 4 double casters
- Central brake acting on 4 wheels
- Power supply-independent battery box with integrated emergency stop button
- Foot keypad for moving the seat up & down
- Extremely easy to move
- Space-saving



EQUIPMENT




Completely height and lengthwise adjustable arm-rests, 350 x 120 mm.



Optional: one-hand height release for direct and high-precision all-around adjustment.
No. 33100000



Adjustable backrest.



Tiltable seat.



Foot keypad enables quick and easy electromotive height adjustment.



Anatomically formed, ergonomic backrest reduces pressure on intervertebral discs and provides ideal support for the lumbar spine.

Seat system evenly distributes your weight as you make the slightest adjustment in position.



Optional: adaptable foot switch plate.

Order No. 33510014

SPECIFICATIONS

Lifting Column

- Drive 3,000 N IPX4
- Lifting height 200 mm

Chassis

- w x l 470 x 535 mm
- Wheels ø 100 mm
- Central brake acting on all four wheels
- Emergency stop button

Weight 45 kg

Max Load 150 kg

Battery unit

- Battery box with 2 batteries (lead-gel), pollution free,

maintenance-free, cadmium-free

2 x 12 V 7.2 Ah

- Universal battery charger:

Input: 100-240 V

Output: 24 V

Seat/backrest

- Min. sitting height 535 mm

- Max. sitting height 735 mm

- Seat adjusting +6°/-9°

- Back adjusting +4°/-18°

Upholstery colors Tundra skai

Standard color anthracite

F6461197

Subject to technical modifications.



surgiForce

ACCESSORIES



FOOTSWITCH PLATE
ADAPTER
Order No. 33170000

ARM REST PAIR
with one-hand release for
neurosurgery (angular),
retrofit kit
Order No. 33500000

surgiTrend

ACCESSORIES



BACKUP BATTERY
Order No. 30470000

FOOTSWITCH PLATE
ADAPTER
Order No. 31050003

surgiLine

ACCESSORIES



BACKUP BATTERY
Order No. 30470000

FOOTSWITCH PLATE
ADAPTER
Order No. 33510014



ARMREST WITH LINEAR
GUIDE
and one-hand release
retrofit kit
Order No. 33510023

ARMREST WITH LINEAR
GUIDE
and one-hand release
factory fitted
Order No. 33510034



ARM REST PAIR
with one-hand release
factory fitted
Art.-Nr. 33100000
retrofit kit
Order No. 33510033

ARMREST WITH LINEAR
GUIDE
and one-hand release
factory fitted
Art.-Nr. 33510035
retrofit kit
Order No. 33510013



ARM REST PAIR
with one-hand release
factory fitted
Art.-Nr. 33100000
retrofit kit
Order No. 33510033

ARMREST WITH LINEAR
GUIDE
and one-hand release
factory fitted
Art.-Nr. 33510035
retrofit kit
Order No. 33510013