



# *supportLine*

INSTRUMENT TABLE



ENJOY QUALITY AND DESIGN

E

# supportLine

Item No.: 3351 0026

The *supportLine* mobile instrument table is designed to be an ergonomic aid to supplement treatment and surgeon's chairs.

**The practical  
little helpers**



## Technical specifications

### Chassis

- Footprint 680x400 mm
- 4 castors ø 50 mm

### Lifting column

- Height adjustment 730 – 985 mm

### Stainless steel tray

660x540 mm

### Upper frame

715x540 mm

### Weight

17,5 kg

### Max. load

10 kg

Subject to technical modifications.

## The perfect addition for your OR

### Pneumatic spring technology/foot pedal

The integrated pneumatic spring technology enables the instrument tray to be stopped and positioned smoothly at a suitable working height for the surgeon. Infinitely variable height adjustment via foot pedal.



# supportLine plus



Item No.: 3351 0045

The *supportLine*'s maneuverability, optimised stability and resistance to tipping, its stainless steel design for the user's operating area and its hard-wearing, hygienic surfaces are designed specifically for versatile use in the operating room. The *supportLine plus* mobile instrument table is a 20 cm higher version.



## Technical specifications

### Chassis

- Footprint 680x400 mm
- 4 castors ø 50 mm

### Lifting column

- Height adjustment 930 – 1,185 mm

### Stainless steel tray

660x540 mm

### Upper frame

715x540 mm

### Weight

17,5 kg

### Max. load

10 kg

Subject to technical modifications..

## Tray

The removable, high-quality stainless steel tray is heat-resistant, easy to clean and disinfect and offers plenty of space for holding surgical instruments with an area of 66 x 54 cm.



## Push handle/double castors

Its wide push handle and maneuverable double castors make the *supportLine* easy to position.



# Upgrade your OR to perfection.



eyeForce



500 XLE



500 XLE comfort



500 XLS



500 ECO



400 XL



motionLine



motionLine-E



surgiForce



surgiForce light



surgiTrend



surgiLine



stayClean wipes



stayClean



staySoft



stayClean roll



supportLine



supportLine plus



assistTrend



assistPro F



assistPro H



easyPro F



easyPro H



easyPro R



easyLine F



easyLine R



easyLine Rs



assistLine

## CONTACT US



**UFSK**  
INTERNATIONAL

**OSYS**

UFSK-International OSYS GmbH  
Kirchhoffstraße 1  
93055 Regensburg  
Germany

Phone +49 (0) 941 78862 15  
Fax +49 (0) 941 78862 35  
Mail [info@ufsk-osys.com](mailto:info@ufsk-osys.com)  
Web [ufsk-osys.com](http://ufsk-osys.com)

**MADE IN GERMANY**  
Design | Construction | Manufacturing