

0.05% TRYPAN BLUE
IN PREFILLED SYRINGE,
0.7 ML

RS-BLUE

0.05% TRYPAN BLUE

TISSUE STAINING



RS-BLUE is 0.05% trypan blue, 0.7 ml, packaged in syringe, for reversible staining of human tissues during surgery. Product warrants a prompt and effective staining of the capsule during cataract surgery, thus allowing a clear visualization of tissues without any risks for patients.

RS-BLUE

Ingredients	<ul style="list-style-type: none"> ▶ trypan blue ▶ sodium chloride ▶ sodium phosphate monobasic ▶ sodium phosphate dibasic ▶ water for injectables
Concentration	0.05% w/V
Hydrogenionic activity (pH)	7.0 – 7.4 pH unit
Osmolality	255-345 mOsm/kg
Endotoxins	< 0,5 EU/ml
Form	liquid
Volume	0.7 ml
Packaging	syringe
Sterilization	moist heat according to EN 556-1 and EN ISO 17665-1
Storage temperature	+15°C / +32°C
Expiry date	36 months
Intended use	reversible tissue staining during ophthalmic surgery
Class	medical device
GMDN code	P 45180
MD registration number	3610/R
Medical device class	Ila
CE mark	CE 0546

Ref.: RSB 002-00

Contents:

- ▶ 6 RS-BLUE prefilled syringes, 0.7 ml
- ▶ 1 insert
- ▶ 12 labels for product traceability

EAN 13 barcode:

8021118000012

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

Tel. +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com