


20	18	16	14	13	12	11	10	9	8	7.5	7	6.5	6	5.5	5.25	5	4.75	4.5	4.25	4	3.75	3.5	3.25	3	2.75	2.5	2.25	2	1.75	1.5	1.25	1	0.75	0.5	0.25	0.12
----	----	----	----	----	----	----	----	---	---	-----	---	-----	---	-----	------	---	------	-----	------	---	------	-----	------	---	------	-----	------	---	------	-----	------	---	------	-----	------	------

6	5	4.5	4	3.5	3	2.75	2.5	2.25	2	1.75	1.5	1.25	1	0.75	0.5	0.25	0.12	MR	FL	⊕	⊕	⊖	⊖	⊙	⊙	RF	GF
---	---	-----	---	-----	---	------	-----	------	---	------	-----	------	---	------	-----	------	------	----	----	---	---	---	---	---	---	----	----

6	5	4.5	4	3.5	3	2.75	2.5	2.25	2	1.75	1.5	1.25	1	0.75	0.5	0.25	0.12	10	12	6	8	4	5	2	3	0.5	1
---	---	-----	---	-----	---	------	-----	------	---	------	-----	------	---	------	-----	------	------	----	----	---	---	---	---	---	---	-----	---

20	18	16	14	13	12	11	10	9	8	7.5	7	6.5	6	5.5	5.25	5	4.75	4.5	4.25	4	3.75	3.5	3.25	3	2.75	2.5	2.25	2	1.75	1.5	1.25	1	0.75	0.5	0.25	0.12
----	----	----	----	----	----	----	----	---	---	-----	---	-----	---	-----	------	---	------	-----	------	---	------	-----	------	---	------	-----	------	---	------	-----	------	---	------	-----	------	------



476 mm (18.7 in)

578 mm (22.8 in)

# Trial Lens Set Comfort

## The classic

The following modules are already included in the Trial Lens Set Comfort

Description	Article no.	Quantity
tl trial lenses +/- sph (basic)	10023545	80
tl trial lenses +/- sph (extension)	10023547	56
tl trial lenses - cyl (basic) <i>axis position 90°</i>	10023556	20
tl trial lenses - cyl (extension) <i>axis position 90°</i>	10023560	16
tl trial lenses special lenses (basic) <i>2 x cross hair, 1 x occluder, 1 x pinhole disc 1.25</i>	10023565	4
tl trial lenses special lenses (extension) <i>coloured lens red, coloured lens green, frosted lens, plano lens, stenopaic slit 1.25 mm, Maddox cylinder red</i>	10023566	6
<b>Trial Lens Set Comfort*</b>	<b>10019887</b>	<b>182</b>

\* uncoated (*this set of trial lenses is also available as coated version*)

*The illustration on the left gives a schematic overview of the trial lens case (article no. 10026133) with the Trial Lens Set Comfort and the optional modules.*

We recommend the following optional modules for expansion

Description	Article no.	Quantity
tl trial lenses +/- sph <i>each 2 x 16/18/20</i>	10023551	12
tl trial lenses + cyl (basic) <i>axis position 90°</i>	10023558	20
tl trial lenses + cyl (extension) <i>axis position 90°</i>	10023562	16
tl trial lenses prism. (basic) <i>base 0°</i>	10023569	5
tl trial lenses prism. (extension) <i>base 0°</i>	10023571	5

**Everything perfectly organized!**

You will find the ideal storage solution for your new trial lenses on **page 14**.