


OCULUS Trial Frames

For your Professional Refraction



OCULUS UB 6

Refraction in perfection

Highlights

- Six slots for each side: four in front, two behind
- Gradation: 2.5°
- PD adjustment R/L: 46 - 80 mm
- Weight: 62 g
- Parallax-free reading of the vertex distance (VD)
- Flexible double earpieces for perfect fit
- Anatomically designed nose pad for more comfort
- Height- and depth-adjustable nose pad
- Delivery content includes three additional nose pads
- Improved thumb wheel for inclination adjustment
- Adjustable arm length
- Made of light metal/plastic



42600	Universal Trial Frame UB 6
42680	Polarizing filters linear (45°/135°) for UB 6 pivotal and reversible (V-/A-position), large viewing area for distance and near vision, full rim, easy fitting with mounting clips
42681	Polarizing filters circular for UB 6 pivotal, large viewing area for distance and near vision, full rim, easy fitting with mounting clips




OCULUS UB 4


Comfort and proven durability

Highlights

- Five slots for each side: three in front, two behind
- Gradation: 5°
- PD adjustment R/L: 46 - 78 mm
- Weight: 64 g
- Parallax-free reading of the vertex distance (VD)
- Flexible bendable earpieces
- Anatomically shaped nose pad
- Height- and depth-adjustable nose pad
- Delivery content includes two additional nose pads
- Thumb wheel for inclination adjustment
- Adjustable arm length
- Made of light metal/plastic



42500	Universal Trial Frame UB 4
45886	Polarizing filters linear (45°/135°) for UB 4 pivotable and reversible (V-/A-position), large viewing area for distance and near vision, full rim, for lateral mount
45888	Polarizing filters linear (90°/180°) for UB 4 pivotable and reversible, click into place for height adjustment for distance or near examination, semi-rimless, for lateral mount
45883	Polarizing filters circular for UB 4 pivotable, large viewing area for distance and near vision, full rim, for lateral mount



OCULUS UB 3+

A tried and tested tradition improved

Highlights

- Five slots for each side: three in front, two behind
- Gradation: 5°
- PD adjustment R/L: 50 - 80 mm
- Weight: 52 g
- Parallax-free reading of the vertex distance (VD)
- Anatomically shaped earpieces
- Delivery content includes one additional nose pad
- Height- and depth-adjustable nose pad
- Thumb wheel for inclination adjustment
- Adjustable arm length
- Made of light metal/plastic



42100	Universal Trial Frame UB 3+
45851	Slip-on pair of rotary linear polarizing filters (45°/135°) by Hegener for UB 3+, reversible

42100 + 45851


OCULUS UB 5

The Universal Trial Frame for small heads


Highlights

- Five slots for each side: three in front, two behind
- Gradation: 5°
- PD adjustment R/L: 40 - 72 mm
- Weight: 48 g
- Parallax-free reading of the vertex distance (VD)
- Flexible bendable earpieces
- Anatomically shaped nose pad
- Height- and depth-adjustable nose pad
- Delivery content includes two additional nose pads
- Thumb wheel for inclination adjustment
- Adjustable arm length
- Made of light metal/plastic
- Delivery content includes five additional coloured knobs

The OCULUS UB 5 is the trial frame for trial lenses of the **BK S Series (tls = trial lens small)** with a rim diameter of **28 mm**.



42550	Universal Trial Frame UB 5
45887	Polarizing filters linear (45°/135°) for UB 5 pivotable and rotatable (V-/A position), large viewing area for distance and near vision, full rim, for lateral mount
45885	Polarizing filters linear (45°/135° and 90°/180°) for UB 5 pivotable and rotatable, large viewing area for distance and near vision, semi-rimless, without midline partition, for lateral mount



Half-eye Trial Frames

Approved Half-eye Trial Frames for adults and children

Highlights

- Three slots for each side
- Gradation: 5°
- PD adjustment R/L:
59 - 67 mm (adults)
54 - 58 mm (children)
- Anatomically shaped earpieces
- Adjustable arm length
- Made of light metal/plastic



Half-eye Trial Frame
rigid nose pad

Half-eye Trial Frame
adjustable nose pad



**Half-eye Trial Frame
for Children**
rigid nose pad

**Half-eye Trial Frame
for Children**
adjustable nose pad

Half-eye Trial Frames	
42580	Half-eye Trial Frame , rigid nose pad
42588	Half-eye Trial Frame , adjustable nose pad
42480	Half-eye Trial Frame for Children , rigid nose pad
42488	Half-eye Trial Frame for Children , adjustable nose pad

Refraction Accessories

Useful accessories for your professional refraction



Pair of polarized filters circular (tl)



Jackson cross cylinder




Demonstrator with Landolt ring



Pair of polarized filters circular (tls)



Refraction bar



Two retinoscopy bars with slider (by Trousseau)

34085	Pair of polarized filters circular (tl), 38 mm
34095	Pair of polarized filters circular (tls), 28 mm

36710	Jackson cross cylinder, +/- cyl 0.12 D, coated
36720	Jackson cross cylinder, +/- cyl 0.25 D, coated
36740	Jackson cross cylinder, +/- cyl 0.50 D, coated
36780	Jackson cross cylinder, +/- cyl 1.00 D, coated
36600	Refraction bar, +/- sph 0.25 und 0.50 D
36605	Refraction bar, +/- sph 0.25 und 0.50 D, coated

46970	Demonstrator with Landolt ring
46980	Demonstrator with E-Hook
48110	Two retinoscopy bars with slider (by Trousseau)

OCULUS Trial Frames

Technical Data

	Weight	PD adjustment range	Arm length
UB 6	62 g (2.19 oz)	46 - 80 mm	92 - 130 mm (3.6 x 5.1 in)
Polarizing filters linear (45°/135°)	12 g (0.42 oz)		
Polarizing filters circular	11 g (0.39 oz)		
UB 4	64 g (2.26 oz)	46 - 78 mm	88 - 124 mm (3.5 x 4.9 in)
Polarizing filters linear (45°/135°)	14 g (0.49 oz)		
Polarizing filters linear (90°/180°)	14 g (0.49 oz)		
Polarizing filters circular	14 g (0.49 oz)		
UB 3+	52 g (1.83 oz)	50 - 80 mm	96 - 132 mm (3.8 x 5.2 in)
Slip-on pair of rotary linear polarizing filters (45°/135°) by Hegener for UB 3+	14 g (0.49 oz)		
UB 5	48 g (1.69 oz)	40 - 72 mm	85 - 121 mm (3.3 x 4.8 in)
Polarizing filters linear (45°/135°)	12 g (0.42 oz)		
Polarizing filters linear (45°/135° or 90°/180°)	15 g (0.53 oz)		
Half-eye Trial Frame with rigid nose pad	25 g (0.88 oz)	59 - 67 mm	94 - 129 mm (3.7 x 5.1 in)
Half-eye Trial Frame with adjustable nose pad	29 g (1.02 oz)	59 - 67 mm	94 - 129 mm (3.7 x 5.1 in)
Half-eye Trial Frame for children with rigid nose pad	25 g (0.88 oz)	54 - 58 mm	94 - 129 mm (3.7 x 5.1 in)
Half-eye Trial Frame for children with adjustable nose pad	29 g (1.02 oz)	54 - 58 mm	94 - 129 mm (3.7 x 5.1 in)

All OCULUS Universal Trial Frames comply with DIN EN ISO 12867.



WWW.OCULUS.DE




The OCULUS QM system is certified in accordance with ISO 13485 (MDSAP) and (EU) 2017/745 (MDR)

OCULUS Optikgeräte GmbH
 Postfach • 35549 Wetzlar • GERMANY
 Tel. +49 641 2005-0 • Fax +49 641 2005-295
 Email: export@oculus.de • www.oculus.de


Find your local OCULUS representative on our website.

Important note

The OCULUS Universal Trial Frames UB 6, UB 4, UB 3+ and all Half-eye Trial Frames are suitable for trial lenses conforming to **DIN EN ISO 9801** with a rim diameter of 38 mm: OCULUS „tl“ trial lens



The OCULUS Universal Trial Frame UB 5 is suitable for trial lenses conforming to **DIN EN ISO 9801** with a rim diameter of 28 mm: OCULUS „tls“ trial lens small



The availability of products and features may vary by country. OCULUS reserves the right to change product specifications and design. All information is valid at the time of printing (02/24).

OC/1895/WZ/EN
 P/SP/028/EN