

ORTOPAD®

ORTOPAD® eye patches are a safe and gentle way to cover the eye and they will remain in place even during energetic play. They are anatomically designed with age-appropriate sizes to ensure they are an excellent fit, preventing adverse light stimuli from reaching the patched eye.




- skin-friendly non-woven bamboo fabric
- breathable
- thin, soft light excluding pad
- ergonomic shape with nasal slits
- hypoallergenic adhesive
- wide selection of motifs


Illustration of the ORTOPAD® outlines in original size:



up to age 2 years



from age 2 to 4 years



upwards of age 4 years



Sample motifs

ORTOPAD[®] SOFT boys

Unique texture eye patches appealing to boys with a soft fleece, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: IN73525
Medium	76 mm x 54 mm	Qty: 50	Stock code: IN73226
Regular	85 mm x 59 mm	Qty: 50	Stock code: IN73227



Sample motifs

ORTOPAD[®] SOFT girls

Unique texture eye patches appealing to girls with a soft fleece, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: IN73225
Medium	76 mm x 54 mm	Qty: 50	Stock code: IN73226
Regular	85 mm x 59 mm	Qty: 50	Stock code: IN73227



Sample motifs

ORTOPAD[®] funpack

Occlusion eye patches with large scale designs for boys and girls. Ergonomic shape, nasal slits and a very skin-friendly adhesive, a lightexcluding centre pad on a soft layer ensure optimal comfort.

Junior	67 mm x 50 mm	Qty: 50	Stock code: US70237
Medium	76 mm x 54 mm	Qty: 50	Stock code: US70238
Regular	85 mm x 59 mm	Qty: 50	Stock code: US70239





Glow in the dark

ORTOPAD
girls

The large scale designs on these occlusion eye patches were created especially to appeal to girls. Five designs with a soft layer, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: IN73225
Medium	76 mm x 54 mm	Qty: 50	Stock code: IN73226
Regular	85 mm x 59 mm	Qty: 50	Stock code: IN73227



ORTOPAD
boys

The large scale designs on these occlusion eye patches were created especially to appeal to boys. Five designs with a soft layer, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: IN73325
Medium	76 mm x 54 mm	Qty: 50	Stock code: IN73326
Regular	85 mm x 59 mm	Qty: 50	Stock code: IN73327



Every box includes colourful stickers for decorating the occlusion eye patches.



ORTOPAD[®]

White occlusion eye patches for use in occlusion treatment. ORTOPAD[®] is made of exceptionally comfortable and breathable fleece, is remarkably soft and especially suitable for use on all sensitive skin types. With light-excluding centre pad and very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: US70021
Medium	76 mm x 54 mm	Qty: 50	Stock code: US70022
Regular	85 mm x 59 mm	Qty: 50	Stock code: US70024



Every box includes colourful stickers for decorating the occlusion eye patches.



ORTOPAD[®] Skin

This beige-coloured occlusion eye patch is extremely gentle to skin. The soft fleece offers outstanding comfort and breathability and comes with a special light-excluding centre pad and very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: US70121
Regular	85 mm x 59 mm	Qty: 50	Stock code: US70124



Information on the adhesive of ORTOPAD® Occlusion Eye Patches

PROPERTIES

The adhesive dressings of the MASTER-AID® product family, to whom ORTOPAD® Eye patches belong, are produced using an adhesive without solvent and latex free. These UV-crosslinked pressure-sensitive adhesives are applied by transfer coating and ensure excellent adhesion to the substrate. For the consumer it is especially important that the product, in comparison to other adhesives, therefore does not cause skin irritation and allergies.

APPLICATION

The pressure-sensitive adhesive is activated only by slight pressure or friction during application. It is then the body temperature which keeps the adhesive properties active. The pressure-sensitive adhesives have optimum adhesive properties, can be removed again painlessly and leave no residue on the skin.

PROTECTION OF THE ENVIRONMENT

The adhesives in use with our products ensure a clean production process. The coating does not release any solvents that are harmful for the environment or the production staff. These features ensure production of unique dressings that excel in environmental protection as well.



ORTOPAD®
eye patches
For best compliance
in daily therapy!

ORTOPAD®

Bamboo as carrier material



FOR THE SAKE OF THE ENVIRONMENT

Working responsibly with nature and its raw materials is very important to us, even more so as we are a supplier of single-use medical devices. We therefore focus on sustainable, resource-saving production when choosing the materials to be used. Therefore our ORTOPAD® eye patches are made from a soft, sustainable bamboo fleece.

Bamboo scores high with its positive eco-balance. The processed raw material simultaneously has very good skin compatibility. The patches are soft, yet durable and allow the skin to breathe optimally.

With a growth rate of up to one metre per day, bamboo is one of the fastest-growing and most renewable raw materials in the world. The special thing about it is that after the harvesting of the stalks the plant remains live, because the rootstock repeatedly sprouts anew. This root system with its widely branching runners also prevents soil erosion.

Compared to traditional woods, bamboo produces much more biomass, thanks to its rapid growth with comparatively low energy input. Bamboo also contributes naturally to reduction of the greenhouse effect. Thanks to its rapid growth, bamboo absorbs a considerable amount of CO₂, whilst producing large amounts of oxygen.

The versatile applicability ultimately make it an ecologically valuable alternative to a variety of resource-scarce raw materials. Therefore we have selected bamboo as the fleece for our eye patches.



The ORTOPAD selection

ORTOPAD
boys

Boxes of 50
 large scale designs




ORTOPAD[®]

Boxes of 50



ORTOPAD[®] Skin

Boxes of 50



ORTOPAD[®] funpack



Boxes of 50

ORTOPAD
girls

Boxes of 50
 large scale designs



Glow in the dark

ORTOPAD[®] SOFT
boys



Boxes of 50

ORTOPAD[®] SOFT
girls



Boxes of 50

ORTOPAD[®] BOYS



Boxes of 50

ORTOPAD[®] GIRLS



Boxes of 50